



NEUBIAS Academy
Online Course

ImageJ/Fiji Macro Language

Anna Klemm - anna.klemm@it.uu.se
BioImage Informatics Facility, SciLifeLab, Sweden

SciLifeLab

With material from:

Fabrice Cordelières, Bordeaux Imaging Center, Bordeaux, France

Giovanni Cardone, MPI Biochemistry, Martinsried, Germany

Ofra Golani, Weizmann Institute of Science, Israel

NEUBIAS material

Outline of This Session



- Biological Data Set and Image Analysis Problem
- How can we "talk" to Fiji? - Macro Recorder, Built-in Macro Function
- Step-by-Step Workflow

The Cell Atlas/Human Protein Atlas



Cell Atlas Aim:


Determine the subcellular location of all cellular proteins.

Experimental Methods:

- Antibody generation against 12.000 human proteins
- Immunostaining, 22 cell lines
- Automated confocal microscopy
- → 82.152 images

Image Analysis Aim:

- Mapping 12.000 human proteins to 30 subcellular structures




Adapted from Thul, P.J. et al. (2017). A subcellular map of the human proteome. Science 356.

The Aim: Quantify Signal Accumulation within the Nuclear Membrane




Image source: Human Protein Atlas

v19.proteinatlas.org/ENSG00000113368-LMNB1




Dataset:

- Subset of The Cell Atlas (Human Protein Atlas)
- 3 color stack: microtubules (magenta), protein detected by antibody (green), nuclei (blue)



*Mean Intensity
Nucleus*



*Mean Intensity
NucMembrane*


$$\text{Ratio} = \frac{\text{Mean Intensity NucMembrane}}{\text{Mean Intensity Nucleus}}$$

The Aim: Quantify Signal Accumulation within the Nuclear Membrane




Image source: Human Protein Atlas

v19.proteinatlas.org/ENSG00000113368-LMNB1



*Mean Intensity
Nucleus*



*Mean Intensity
NucMembrane*

$$\text{Ratio} = \frac{\text{Mean Intensity NucMembrane}}{\text{Mean Intensity Nucleus}}$$

Results		
File	Edit	Font Results
	Label	Mean
1	signal:0001-0041	26957.406
2	signal:0002-0291	20013.618
3	signal:0003-0320	38092.890
4	signal:0004-0670	18741.716
5	signal:0005-0696	19940.679
6	signal:0006-0677	16445.010
7	signal:0007-1198	20677.366
8	signal:0008-1168	20005.914
9	signal:0009-1250	24444.675
10	signal:0010-1457	20037.129
11	signal:0011-1651	26454.839
12	signal:0012-1788	5380.207
13	signal:0013-1832	24655.655
14	signal:0014-1945	15790.756
15	signal:0015-1935	18773.019
16	signal:0016-2019	17444.867
17	signal:0017-2021	9947.993

Exercise: Find a Workflow




Task: Create selections (ROIs) around nuclei and the nuclear membrane

- Create selections for each nucleus: Threshold, Connected Component Analysis
- Try for one nucleus: create a selection around the nuclear membrane.

*Image source: Human Protein Atlas
v19.proteinatlas.org/ENSG00000113368-LMNB1*

Image:

- CellAtlas_Subset/711_D6_1.tif





What would be your workflow?

Create a selection for each nucleus:

- 1) Split channels
- 2) Prefilter and segment nuclei
- 3) Run „Analyze Particles“ to get ROIs of single Nuclei

Try for one nucleus: ROI of the nuclear membrane

- 1) Select the ROI of one nucleus
- 2) Decrease the size of the ROI by
Edit>Selection>Enlarge
- 3) *Edit>Selection>Make Band...* to create ROI of the envelope





We now want to automatize this process



IMAGEJ MACRO LANGUAGE

Recorder




Open record window:

- Plugins > Macros > Record...
- Choose language – if needed




Let's try with only one command: Split Channels



Recorder

Recorder:

- Discover commands
- Window can be edited, copied, pasted, cut etc.





Script Editor

Script Editor:

- For building a script
- Has color-coding, code-completion, run-option etc.

How to start a new Script



Importance of recorder



- Discover commands
- Get arguments for specific functions

But also:

- Record and save your workflow – for documentation and reproducibility!

Exercise: Record the Workflow



Task:

- Open the recorder and record the discussed workflow.
- Discuss open questions within the group.

Image:

- Any image from CellAtlas_Subset

The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name

Filter nuclear image and make binary image

Retrieve nuclei boundaries

Retrieve the nuclear envelope boundary

Measure

The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name



Programming Basics I

VARIABLES

Variables: definition



```
*New_.ijm
1 totalArea = 100;
2 fileName = "wildtype.tif";
3 description = "Launching the script...";
4 thereAreCells = true;
```

- Can hold numbers or phrases/strings, but only one at a time
- Used whenever a value is used many times inside the script
- You define a variable by assigning it some content
- Variable name is on the left followed by an equal sign followed by the item (or items) being assigned
- Variable names can only start with characters

Numeric Variables: assignment




```
1 totalArea = 0;
2 numPixelsCell1 = 154;
3 numPixelsCell2 = 190;
4 pixelSize = 0.350;
5 totalArea = (numPixelsCell1 + numPixelsCell2)*pixelSize;
6 print(totalArea);
```

Run Kill Show Errors Clear

Started New_.ijm at Fri Mar 17 13:24:57 CET 2017

- Content of numeric variables can be modified using mathematical operations
- After an assignment, the previous content (if any) is forgotten
- Good practice is to use **d2s** (decimal to string) when printing numbers
e.g. `print(d2s(totalArea));`



String Variables: concatenation



```
string_concatenation.ijm
1 text1 = "a";
2 text2 = "Hello";
3 text3 = "Hello everybody!";
4 text4 = " ";
5
6 text5 = text1 + text2;
7 print(text5);
8 print(text1 + text2);
9 print(text2+ "world!");
10 print(text2 + text4 + "world!");
```

```
string_numbers_concatenation.ijm
1 number1 = "2";
2 number2 = "3";
3 print(number1+number2);
4
5 number3 = 2;
6 number4 = 3;
7 print(number3+number4);
8
9 text = "image";
10 print(text + number3);
```

```
*New_.ijm
File Edit Language Templates Run Tools Tabs
*New_.ijm
1 fileName = "wt007.tif";
2 print(fileName);
Run Kill Show Errors Clear
```

What about this?



string_concatenation.ijm



string_numbers_concatenation.ijm

The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name

Introduction of the next Exercise



- Define a string variable with the name of the image
- Build selectWindow using the variables

```
run("Split Channels");  
  
selectWindow("C1-711_D6_1.tif");  
rename("nuclei");  
  
selectWindow("C2-711_D6_1.tif");  
rename("signal");
```

Exercise: Usage of Variables



- Replace the **highlighted text** using the variable *title*.
- Check out slide 19 (String Variables: concatenation) for help.

```
title = "711_D6_1.tif";  
  
run("Split Channels");  
selectWindow("C1-711_D6_1.tif");  
rename("nuclei");  
  
selectWindow("C2-711_D6_1.tif");  
rename("signal");
```


Structuring the Code



Technical point:

Structuring the code using comments

- Comments are non-interpreted elements of code
- They help structure the code
- They help collaborators interpret the original analyst's intentions
- Comments are introduced either by // or surrounded by /* */:

```
//This is a short comment
```

```
/*
```

```
    This is a very long comment, spanning over  
    multiple lines, allowing line breaks
```


```
*/
```



We now know what variables are but..

HOW TO GET THE NAME OF AN IMAGE AUTOMATICALLY?

Built-in Macro Functions



- Not everything is recorded. Much more functions can be found at:
<https://imagej.nih.gov/ij/developer/macro/functions.html>
- Tip: do a page-search (CTRL+F)

Step 1:

Normalise the data name



Solution

```
Step_01_SplitAndRename.ijm
File Edit Language Templates Run Tools Tabs Options
Step_05_Measure.ijm Step_05_Measure.ijm Step_01_SplitAndRename.ijm
1
2
3
4
5
6
7
8 //Step1: Getting image information + Normalise the data name
9 //get general information
10 title = getTitle();
11
12
13 //split channels and rename them
14 run("Split Channels");
15 selectWindow("C1-" + title);
16 rename("nuclei");
17 selectWindow("C2-" + title);
18 rename("signal");
19
20
21
22
```



Catch-up/[Step_01_SplitAndRename.ijm](#)

More about the Built-in Macro Functions



Stack.setChannel(1);

Function with input

getTitle();

Function with output;

nameOfMyImage = getTitle(); output is assigned to a variable

getDimensions(width, height, channels, slices, frames);

Output is assigned to variables within the brackets

Exercise: Built-in Macro Functions



Task 1: Catch-up with the script

- include the getTitle() function.





Task 2: explore the built-in macro functions.

- open a new script, set language to IJ1 macro
- What happens when you run getDimensions(channels, height, width, slices, frames)?
- Use the print() function to explore the content of the variables channels and width.
- Check the usage of the getDimensions function either using code autocompletion or on the “built-in macro function” website.

The script editor supports you with colors and error messages.




Color-coding in Script-Editor




Read the error messages!

Auto-completion helps to avoid common mistakes



Workflow: How will we tackle the problem ?



The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Filter nuclear image and
make binary image

Retrieve nuclei boundaries

Step 2 + 3: *Pre-process the nuclear image*



```
Step_02_03_Preprocess_AnalyzeParticles.ijm
File Edit Language Templates Run Tools Tabs Options

Step_02_03_Preprocess_AnalyzeParticles.ijm

21 //Step2: Prefilter nuclear image and make binary image
22 selectWindow("nuclei");
23 //preprocessing of the grayscale image
24 run("Median...", "radius=8");
25 //thresholding
26 setAutoThreshold("Huang dark");
27 setOption("BlackBackground", true);
28 run("Convert to Mask");
29 //postprocessing of binary image
30 run("Fill Holes");
31
32
33 //Step3: Retrieve the nuclei's boundaries
34 num = getNumber("minimum size", 2000 );
35 selectWindow("nuclei");
36 run("Analyze Particles...", "size=" + num + "-Infinity add"); //add
37
```



Step_02_03_Preprocess_AnalyzeParticles.ijm

Exercise: Preprocessing, Asking for User-Input




Task 1: Catch-up with the script

As shown in the slide before:

- Insert the preprocessing and segmentation steps to your code (median filter, thresholding, fill holes)
- Use the `getNumber` function to ask the user for a minimum size of the nuclei in pixel

Workflow: How will we tackle the problem ?



The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name

Filter nuclear image and make binary image

Retrieve nuclei boundaries

Step 4:

Retrieve the nuclear envelope's boundaries



```
Recorder
Record: Macro Name:
roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");
```

Selects only one specific ROI

Shrinks the ROI

Creates donut-Shape ROI

Replaces the active ROI with the donut

We need to loop this code over all ROIs!



Programming Basics

CONTROL STRUCTURES FOR-LOOPS

Programming Basics: Control Structures



Technical point

Loops

Definite loop

- Known number of iterations
- 3 arguments:
 - Initialisation
 - Condition for loop entry, as a boolean
 - Iteration

```
for(i=0; i<10; i++){  
    //Instruction 1  
  
    //Instruction 2  
  
    //Instruction 3  
}
```

Indefinite loops

A priori

- Test performed **BEFORE** instructions are executed

```
i=0;  
while(i<10){  
    //Instruction 1  
    //Instruction 2  
    //Instruction 3  
  
    i++;  
}
```

A posteriori

- Instructions always executed at least once
- Test performed **AFTER** instructions have been executed

```
i=0;  
do{  
    //Instruction 1  
    //Instruction 2  
    //Instruction 3  
  
    i++;  
} while(i<10)
```

Programming Basics: For-Loops



```
ForLoop_Example1.ijm
File Edit Language Templates Run Tools Tabs
ForLoop_Example1.ijm
1 //use i only as condition for the for-loop
2 for (i=0; i<10; i++){
3     print("Neubias is great!");
4 }
5
6 //use i additionally for calculations
7 for (i=0; i<10; i++){
8     result = i * 10;
9     print(result);
10 }
11
Run Kill Show Errors Clear
```



Catch-Up/ForLoop_Example1.ijm

```
*ForLoop_Example2.ijm (Running)
File Edit Language Templates Run Tools Tabs
ForLoop_Example1.ijm *New_.ijm *ForLoop_Example2.ijm
1 //different writing for i++
2 for (i=0; i<10; i+=1){
3     print(i);
4 }
5
6 //use a different increment
7 for (i=0; i<10; i+=2){
8     print(i);
9 }
10
11
Run Kill Show Errors Clear
```



Catch-Up/ForLoop_Example2.ijm

For-Loops Exercise



Find **four different ways** to modify the code below to print „Hello!“ **only once** instead of 10 times.

```
for (i=0; i<10; i++){  
    print("Hello!");  
}
```

Exercise: Looping



Task:


Modify the recorded code to apply it
to all ROIs of the ROI-Manager.

```
roiManager("Select", 7);  
run("Enlarge...", "enlarge=-4");  
run("Make Band...", "band=7");  
roiManager("Update");
```

Hints:

- Search for an built-in function giving the **total number of ROIs**.
- Index of the first ROI is **0**!
- **Close the ROI-manager** each time before re-running your code again and again.

Close the ROI-manager before Running again



For-Loop run over all ROIs, also “old” ones
→ Close the ROI-manager!

We will **automatically** close the ROI-manager soon...

Step 4: *Retrieve the nuclear envelope's boundaries*



Solution


```
*Step_04_ForLoop.ijm
File Edit Language Templates Run Tools Tabs Options

*Step_04_ForLoop.ijm
40
41
42 //Step4: Retrieve the nuclear envelope's boundaries and save them in the ROI-Manager
43 numberOfNuclei = roiManager("count");
44 for(i=0; i<numberOfNuclei; i++){
45     roiManager("Select", i);
46     run("Enlarge...", "enlarge=-4");
47     run("Make Band...", "band=7");
48     roiManager("Update"); //original nucleus-ROI is replaced by nuclear envelope ROI
49 }
50
51
```



Catch-up/Step_04_ForLoop.ijm

Workflow: How will we tackle the problem ?



The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name

Filter nuclear image and make binary image

Retrieve nuclei boundaries

Retrieve the nuclear envelope boundary

Measure







Step 5:

Make the measurement and save the result

Hint:

If no ROI is selected in the ROI Manager, `roiManager("Measure");` will measure **all ROIs**.





Empty ROI Manager and Empty Results window

Make sure to start with an empty ROI-manager and an empty Results table!

```
roiManager("reset");  
run("Clear Results");
```

```
Step_06_CleanStart.ijm  
File Edit Language Templates Run Tools Tabs Options  
Recorded.ijm Step_06_CleanStart.ijm Step_07_batchProcessing.ijm  
1 /*  
2  * NEUBIAS Academy  
3  * ImageJ/Fiji Macro Language  
4  * anna.klemm@it.uu.se - BioImage Informatics Facility @ScilifeLab  
5  * April 2020  
6  */  
7  
8 //clean-up to prepare for analysis  
9 roiManager("reset");  
10 run("Clear Results");  
11  
12 //Step1: Getting image information + Normalise the data name  
13 //get general information  
14 title = getTitle();  
15  
16
```

Exercise: Catch-up



Task:

- Insert the missing steps to your code
- For saving: Insert a path fitting to your local computer.

```
Step_05_Measure.ijm
File Edit Language Templates Run Tools Tabs Options
Step_02_03_Preprocess_AnalyzeParticles.ijm Recorded.ijm Step_05_Measure.ijm
48 //Step 5: Measure signal in nuclear envelope's boundaries and save the result
49 run("Set Measurements...", "mean display redirect=None decimal=3"); //define the
50 selectWindow("signal");
51 roiManager("deselect"); //ensures that no ROI is selected
52 roiManager("Measure"); //measures active ROI or - if no ROI is selected - all RO
53 // Save results
54 saveAs("results", "C:/Users/Anna/Desktop/Neubias_output/Results.csv");
55
```

Attention: Use / or \




Catch-Up/Step_05_Measure.ijm



EXTRA-STEPS

Extra I: Batch Processing



Extra I: Batch processing – getFileList



- Example script:

fileList is
an **array**

```
input_path = getDirectory("input files");  
fileList = getFileList(input_path);  
  
for (i=0; i<fileList.length; i++){  
    print(fileList[i]);  
}
```


number of images
in the array
(=directory)



Extra_02_getFileList_Example.ijm

Prints filename on position i in fileList.
Attention: 1st position has index 0!

Arrays are multi-compartment boxes with a label



- Can hold numbers or phrases/strings, different type for each compartment/cell
- Can only be 1D
- Used to keep similar data together (e.g. points in plot, measures at different slices, etc.)
- Can be defined everywhere inside a script
- Index starts at 0 and ends at N-1, where N is the size of the array

Arrays: definition



Two ways to initialize an array

```
fruit = newArray(4);  
fruit[0] = "apple";  
fruit[1] = "orange";  
fruit[2] = "apricot";  
fruit[3] = "pear";
```

```
fruit = newArray("apple", "orange",  
                "apricot", "pear");
```

newArray(size) returns a new array containing *size* elements, all with value 0

How to access and print the content of an array

```
for (i =0; i<fruit.length;  
i++){  
    print(fruit[i]);  
}
```

```
Array.print(fruit);
```

fruit.length

returns the number of elements composing the array

Array.print

is one of the functions available for arrays: see the documentation for a list of them

Break-out Exercise: Batch-processing



Task:

- Use the template below to loop your code over all files in the input folder.

```
input_path = getDirectory("input files");  
fileList = getFileList(input_path);
```

```
for (f=0; f<fileList.length; f++){
```

```
    //Clean-up to prepare for next image
```

```
    roiManager("reset");
```

```
    run("Close All");
```

```
    run("Clear Results");
```

```
    open(input_path + fileList[f]);
```

```
    print(input_path + fileList[f]); //displays file that is processed
```

```
    // rest of the code
```

```
    // (...)
```

```
    saveAs("results", "C:/Users/Anna/Desktop/"+title+"_results.xls");
```

```
}
```



Step_06_batchProcessing.ijm





Programming Basics

CONTROL STRUCTURES IF-STATEMENTS

Extra II: conditional statements – IF/ELSE IF/ELSE



- Aim:
 - To ask for the “minimum size” for *Analyze Particles..* only for the first image of a set.

```
//Step3: Retrieve the nuclei's boundaries  
num = getNumber("minimum size", 2000);  
run("Analyze Particles...", "size=" +num+ "-Infinity add");
```

Programming Basics: Conditional statements



Technical point *Conditional execution*

Check one condition

- Only performs the operation if a condition is fulfilled
- The condition might result from a test:
 - Regular: <, >, <=, >=
 - Equality: ==
 - Difference: !=
- The result of the test could be stored in a variable (boolean)

```
i = 0;
if(i > 10){
    //Instruction 1
    //Instruction 2
    //Instruction 3
}
```

...one alternative

- Branch on 2 sets of instructions

```
i = 0;
if(i > 10){
    //Instruction 1
    //Instruction 2
}else{
    //Instruction 3
    //Instruction 4
}
```

...multiple alternatives

- Branch or more than 2 alternatives

```
i = 0;
if(i > 10){
    //Instruction 1
}else if(i < 8) //Instruction 2
}else if(i < 6 ) //Instruction 3
}else //Instruction 4
```

Programming Basics: If-statements



```
1 area = 5; //area in  $\mu\text{m}^2$ 
2
3 if (area < 20){
4     print("Not a healthy Hela-cell!");
5 }
```

```
1 //if-statement inside a loop
2 for (i = 0; i < 10; i++) {
3     print(i);
4     if(i==4){
5         print("Half way done!");
6     } else if(i==9){
7         print("Done!");
8     }
9 }
```



IfStatement_Example.ijm

Break-out Exercise: If-statement



Task:

To ask for the “minimum size” for *Analyze Particles..* only for the first image of a set.

```
//Step3: Retrieve the nuclei's boundaries  
num = getNumber("minimum size", 1000);  
run("Analyze Particles...", "size=" +num+ "-Infinity add");
```


Extra II: conditional statements – IF/ELSE



```
Step_07_batchProcessing_if.ijm
File Edit Language Templates Run Tools Tabs Options
Step_07_batchProcessing_if.ijm
7 //batch processing
8 input_path = getDirectory("Choose image folder");
9 fileList = getFileList(input_path);
10
11 for (f=0; f<fileList.length; f++){ //loops over all images in the giv
12     open(input_path + fileList[f]);
13     print(input_path + fileList[f]); //displays file that is processe
14
15     if(f==0){
16         num = getNumber("minimum size", 1000); //minimum particle siz
17     }
18     //Step1: Getting image information + Normalise the data name
19     //get general information
20     title = getTitle();
21
22     //remove scale to work with pixels
23     run("Set Scale...", "distance=0 known=0 pixel=1 unit=pixel");
24
25     //split channels and rename them
26     run("Split Channels");
```



Step_07_batchProcessing_if.ijm

The recorded workflow



Cleaned-up recorder:



Catch-up/Recorded.ijm

```
Recorder
Record: Macro Name: Macro.ijm Create ?
run("Split Channels");
selectWindow("C1-711_D6_1.tif");
rename("nuclei");
selectWindow("C2-711_D6_1.tif");
rename("signal");

selectWindow("nuclei");
run("Median...", "radius=8");
setAutoThreshold("Huang dark");
setOption("BlackBackground", false);
run("Convert to Mask");
run("Fill Holes");

run("Analyze Particles...", "size=2000-Infinity add");
roiManager("Show None");

roiManager("Select", 7);
run("Enlarge...", "enlarge=-10");
run("Make Band...", "band=10");
roiManager("Update");

run("Set Measurements...", "mean display redirect=None decimal=3");
selectWindow("signal");
roiManager("Select", 7);
run("Measure");
```

Normalise the data name

Filter nuclear image and make binary image

Retrieve nuclei boundaries

Retrieve the nuclear envelope boundary

Measure



Readability:

- Use meaningful variable names
- Assign variable at the top of a script if the variable is used widely, or as close to where it is used as possible
- Comment your code: for you and others

Reproducibility:

- Add Initialization code: close windows, reset roiManager, reset Results table ...
- Save quality control files, e.g. save the ROI manager
- Use file names that refer to the original files
- Save the parameters used with the other results.
- Save the macro itself or document its version
- Consider sharing your macro and parameters as Supplementary Information

Where to continue (I).




There are many more very useful built-in macro functions.

Check out [SneakPreview_calculateRatio.ijm](#) to

- Get all measurements: *Mean Intensity NucMembrane*, *Mean Intensity Nucleus* and *Ratio*
- See how to customize the result window
- Get basic statistics of ROIs using `getStatistics`
- `waitForUser`: easy method to check intermediate results, debug, or wait for user interaction.

Results					
File Edit Font Results					
	image	mean Int nuclei	mean Int nuclear membrane	ratio	minimum size nuclei
1	931_H7_2.tif	34896.065	45857.361	1.314	700
2	931_H7_2.tif	28762.512	37526.908	1.305	700
3	931_H7_2.tif	27779.874	49714.626	1.790	700
4	931_H7_2.tif	31307.901	19444.399	0.621	700

Where to continue (II).



Other resources:

imagej.net/Introduction_into_Macro_Programming

forum.image.sc: Forum thread linked to this video

Chapter 3, ImageJ Macro Language
(free download)

Data resource



Raw images (tif) were provided by The Human Protein Atlas.

<https://www.proteinatlas.org/humanproteome/cell>

Thul, P.J. et al. (2017). A subcellular map of the human proteome. *Science* 356.